

# Instructions for use


## Product description Ultrasonic Cleaner

Date of revision: 10/03/23 v1.5


topdental™  
Bossklein®

Enzyme based concentrate for the pre-cleaning and decontamination of the non-porous surfaces of hand instruments and forceps, rotary instruments, mirrors etc. Designed for use in the Ultrasonic bath. Specially formulated for use in the Medical and Dental environment.


### 5L bottle



1. Add 30ml ± 1ml of the concentrate to 970ml tap water (3% solution).




2. In the ultrasonic bath completely immerse the instruments into the solution ensuring that all surfaces are covered by the solution, add hollow objects in such a way that air escapes.




15 mins

3. Run the ultrasonic cleaner for a minimum of 15 minutes and according to the manufacturers instructions.




4. After immersion, remove the instruments and rinse them thoroughly prior to terminal sterilisation or high level disinfection.




5. Prepare a fresh solution daily or sooner if solution is very dirty. For professional use only. See product labels and MSDS for further information.

### 1L dispensing bottle



1. Add 30ml ± 1ml of the concentrate to 970ml tap water (3% solution) using the dosing device on the bottle.




2. In the ultrasonic bath completely immerse the instruments into the solution ensuring that all surfaces are covered by the solution, add hollow objects in such a way that air escapes.




15 mins

3. Run the ultrasonic cleaner for a minimum of 15 minutes and according to the manufacturers instructions.



4. After immersion, remove the instruments and rinse them thoroughly prior to terminal sterilisation or high level disinfection.



5. Prepare a fresh solution daily or sooner if solution is very dirty. For professional use only. See product labels and MSDS for further information.

**Composition:** Aqua, Propan-2-ol, Didecyltrimethylammonium chloride, Sodium n-lauroylsarcosinate, Subtilisin, Cl 19140.

**Warning.** H226: Flammable liquid and vapour. H315: Causes skin irritation. H319: Causes serious eye irritation. H336: May cause drowsiness or dizziness. P302+350: IF ON SKIN: Gently wash with plenty of soap and water. P305+351+338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P280: Wear protective gloves/protective clothing/eye protection/face protection. P261: Avoid breathing dust/fumes/gas/mist/vapours/spray. P304+340: IF INHALED: Remove person to fresh air and keep comfortable for breathing. **Storage.** P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P403+P233: Store in a well-ventilated place. Keep container tightly closed. P235: Keep cool. **Contraindications.** Incompatible materials: Aluminium, copper. Not recommended for plastics and rubbers. **Disposal.** Do not discharge concentrate into drains or rivers.

**Intended Users.** The products are used within the dental surgery by the dentist or dental assistant as part of their general cleaning/disinfecting routines. Patients are usually not present during use, but all patients that are seen in the surgery contribute to the need for the products.

Any serious incident related to this medical device must be reported to the manufacturer and the competent authorities.

UFI: 0610-EQJV-7002-VXNE



DIRECTA AB  
P.O. Box 723, Finvids väg 8 SE-194 27  
Upplands Väsby, Sweden



For additional languages  
see [www.topdental.co.uk](http://www.topdental.co.uk)



Warning.

DirectaDentalGroup  
DIRECTA topdental ORSING  
parkell KOHLER RÖNVIG  
TrollDental Sendoline CONTACZ  
physics Forcerip ceramir