
ParaPost® X System

ParaPostXP
ParaPostPlus
ParaPost
One-Office-Visit Technique

ParaPostXH
One-Office-Visit Technique

ParaPostXP
ParaPost
Casting Technique
(Direct & Indirect)

Instructions for Use

CE
0344

coltène
whaledent

EN	Instructions For Use	3
DE	Gebrauchsanweisung	12
FR	Mode d'emploi	22
ES	Instrucciones de uso	32
IT	Istruzioni per l'uso	42
NL	Gebruiksaanwijzing	52
SV	Användarinstruktioner	62
DA	Brugsanvisning	72
FI	Käyttöohjeet	81
NO	Bruksanvisninger	90
PT	Instruções para uso	100
EL	Οδηγίες Χρήσης	110

ParaPostXP
ParaPost Plus
ParaPost

One-Office-Visit Technique

For illustrations see page 128

Note that titanium posts have virtually the same radiopacity as composite and gutta-percha. Careful scrutiny of the radiograph will allow for the differentiation of post from gutta percha and composite.

We urge the use of Hygenic[®] Rubber Dam with this procedure.

1. After endodontic therapy has been completed, prepare tooth for restoration as if pulp were intact. The preparation should include at least 1.5 mm of sound tooth structure all around the circumference of the preparation, apical to the core, for desired ferrule effect (**Fig. 01**).* If that amount of tooth structure is not available, it should be created by periodontal crown lengthening and/or orthodontic extrusion.*
2. Use a radiograph to determine the appropriate diameter and depth of post space preparation. Keep in mind that sufficient root wall thickness must be maintained (1 mm minimum)* to prevent perforation or

weakening of root walls, and at least 4-5 mm of gutta-percha must be left intact to protect the apical seal.

3. Remove gutta-percha to preplanned depth with a Gates-Glidden drill, Peeso reamer and/or hot instrument. Radiographic verification is recommended. If the canal has been obturated with a silver cone, remove it and reseal with Hygenic Gutta-Percha before preparing the post space.
4. To begin paralleling the post space, select the ParaPost Drill diameter which corresponds to the last Gates-Glidden drill used in Step 3 (**Fig. 01**). (See “Approximate Diameter Comparisons” chart, **Fig. A**). To help maintain the preplanned depth of the preparation, ParaPost Drills have been pre-marked to identify depths of 7, 9 and 11 mm from the apical end of the drill.

ParaPost Drills can be used either manually with the Universal Hand Driver or with a slow-speed contra-angle. When a contra-angle is used, the drill must be kept in continuous clockwise rotation until it has been completely removed from the tooth. This will minimize the risk of the drill seizing in the post space. Remove all weakened or unsupported tooth

- structure. Remove any tooth debris from canal by irrigating post space with water spray.
5. Sequentially step up to the next larger ParaPost Drill until the preplanned diameter and depth are achieved (**Fig. 02**).
 6. Use a cylindrical diamond or carbide bur to prepare an anti-rotational box (**Fig. 03**).
 7. Select the ParaPost post that corresponds to the last drill used to prepare the post space. Remove head portion from the post. The head portion is for identification purposes only. Insert post into post space (**Fig 04**). Check for occlusal clearance. Remove the post; shorten if necessary.
 8. To prepare for cementation, clean and rinse the canal (**Fig. 05**) and dry with Hygienic Paper Points. Use the dental cement of your choice according to the manufacturer's instructions.
 9. Use a Lentulo spiral to coat the walls of the post space with cement (**Fig. 06**).
 10. Coat the post shaft with cement and slowly insert the post to full depth, allowing excess cement to vent (**Fig. 07**).
 11. If necessary, to obtain additional occlusal clearance after the cement has set, make any adjustments needed to the post head while holding the post with a hemostat (**Fig. 08**).
 12. After the cement has set, remove any excess before completing the core and fabricating the final restoration (**Fig. 09**).

Cleaning and Sterilization

Drills: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush; thoroughly rinse and dry. Place drills, with a sterilization indicator, into an autoclave bag and seal. Sterilize in autoclave at 132°C for 12 minutes. Sterilizing bags should not touch the walls of the autoclave. Once the sterilization monitor indicates that sterilization was achieved, place autoclave bag, with the drills inside, in storage until ready for use. To decontaminate, soak in a 70% isopropyl alcohol solution for 10 minutes.

Posts: Prior to each use, place posts, with a sterilization indicator, into an autoclave bag and seal. Sterilize in an autoclave at 130 - 135 degrees C for 12-15 minutes, following standard procedures.

Universal Hand Driver: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush. Sterilize by autoclave or chemical vapor sterilization.

Warning Notice: DO NOT BEND POSTS.

**References and/or literature available upon request.*

ParaPostXH

One-Office-Visit Technique

For illustrations see page 129

Note that titanium posts have virtually the same radiopacity as composite and gutta-percha. Careful scrutiny of the radiograph will allow for the differentiation of post from gutta percha and composite.

We urge the use of Hygenic[®] Rubber Dam with this procedure.

1. After endodontic therapy has been completed, prepare tooth for restoration as if pulp were intact. The preparation should include at least 1.5 mm of sound tooth structure all around the circumference of the preparation, apical to the core, for desired ferrule effect (**Fig. 10**)* If that amount of tooth structure is not available, it should be created by periodontal crown lengthening and/or orthodontic extrusion.*
2. Use a radiograph to determine the appropriate diameter and depth of post space preparation. Keep in mind that sufficient root wall thickness must be maintained (1 mm minimum)* to prevent perforation or weakening of root walls, and at least 4-5 mm of gutta-percha must be left intact to protect the apical seal.

- Remove gutta-percha to preplanned depth with a Gates-Glidden drill, Peeso reamer and/or hot instrument. Radiographic verification is recommended. If the canal has been obturated with a silver cone, remove it and reseal with Hygenic Gutta-Percha before preparing the post space.
- To begin paralleling the post space, select the ParaPostX Drill diameter which corresponds to the last Gates-Glidden drill used in Step 3 (**Fig. 10**). (See “Approximate Diameter Comparisons” chart, **Fig. A**). To help maintain the preplanned depth of the preparation, ParaPostX Drills have been pre-marked to identify depths of 7, 9 and 11 mm from the apical end of the drill.

ParaPostX Drills can be used either manually with the Universal Hand Driver or with a slow-speed contra-angle. When a contra-angle is used, the drill must be kept in continuous clockwise rotation until it has been completely removed from the tooth. This will minimize the risk of the drill seizing in the post space. Remove all weakened or unsupported tooth structure. Remove any tooth debris from canal by irrigating post space with water spray.
- Sequentially step up to the next larger ParaPostX Drill until the preplanned diameter and depth are achieved (**Fig. 11**).
- Use a cylindrical diamond or carbide bur to prepare an anti-rotational box (**Fig. 12**).
- Select the ParaPostXH post that corresponds to the last drill used to prepare the post space. Insert post into post space (**Fig. 13**). Check for occlusal clearance. Remove post; if post has to be shortened, make adjustments from apical end (**Fig. 14**). Re chamfer the end of the post to its original shape.
- To prepare for cementation, clean and rinse the canal (**Fig. 15**) and dry with Hygenic paper points. Use the dental cement of your choice according to the manufacturer’s instructions.
- Use a Lentulo spiral to coat the walls of the post space with cement (**Fig. 16**).
- Coat the post shaft with cement and slowly insert the post to full depth, allowing excess cement to vent (**Fig. 17**).

11. If necessary, to obtain additional occlusal clearance after the cement has set, make any adjustments needed to the post head while holding the post with a hemostat (**Fig. 18**).
12. After the cement has set, remove any excess before completing the core and fabricating the final restoration (**Fig. 19**).

Cleaning and Sterilization

Drills: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush; thoroughly rinse and dry. Place drills, with a sterilization indicator, into an autoclave bag and seal. Sterilize in autoclave at 132°C for 12 minutes. Sterilizing bags should not touch the walls of the autoclave. Once the sterilization monitor indicates that sterilization was achieved, place autoclave bag, with the drills inside, in storage until ready for use. To decontaminate, soak in a 70% isopropyl alcohol solution for 10 minutes.

Posts: Prior to each use, place posts, with a sterilization indicator, into an autoclave bag and seal. Sterilize in an autoclave at 130 - 135 degrees C for 12-15 minutes, following standard procedures.

Universal Hand Driver: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush. Sterilize by autoclave or chemical vapor sterilization.

Warning Notice: DO NOT BEND POSTS.

**References and/or literature available upon request.*

ParaPostXP
ParaPost

Casting Technique - Direct and Indirect

For illustrations see page 130-131

We urge the use of Hygenic[®] Rubber dam with these procedures.

1. After endodontic therapy has been completed, prepare tooth for restoration as if pulp were intact. The preparation should include at least 1.5 mm of sound tooth structure all around the circumference of the preparation, apical to the core, for desired ferrule effect (**Fig. 20**).^{*} If that amount of tooth structure is not available, it should be created by periodontal crown lengthening and/or orthodontic extrusion.^{*}
2. Use a radiograph to determine the appropriate diameter and depth of post space preparation. Keep in mind that sufficient root wall thickness must be maintained (1 mm minimum)^{*} to prevent perforation or weakening of root walls, and at least 4-5 mm of gutta-percha must be left intact to protect the apical seal.

3. Remove gutta-percha to preplanned depth with a Gates-Glidden drill, Peeso reamer and/or hot instrument. Radiographic verification is recommended. If the canal has been obturated with a sliver cone, remove it and reseal with Hygenic Gutta-Percha before preparing the post space.
4. To begin paralleling the post space, select the ParaPost Drill which corresponds to the last Gates-Glidden drill used in Step 3 (**Fig. 20**). (See “Approximate Diameter Comparisons” chart, **Fig. A**). To help maintain the preplanned depth of the post space, ParaPost Drills have been premarked to identify depths of 7, 9 and 11 mm from the apical end of the drill.

ParaPost Drills can be used either manually with the Universal Hand Driver or with a slow-speed contra-angle. When a contra-angle is used, the drill must be kept in continuous clockwise rotation until it has been completely removed from the tooth. This will minimize the possibility of the drill seizing in the post space. Remove all weakened or unsupported tooth structure. Remove any tooth debris from canal by irrigating post space with water spray.

5. Sequentially step up to the next larger ParaPost Drill until the previously determined diameter and depth are achieved (**Fig. 21**).
6. Use a cylindrical diamond or carbide bur to prepare an anti-rotational box without undercuts (**Fig. 22**).

Direct Only

7. Place into the prepared post space either a ParaPost Burnout Post or a ParaPost No-Ox Precious Alloy Post, that corresponds to the largest ParaPost Drill used to prepare the post space (**Fig. 23**). If necessary, shorten the Burnout Post or No-Ox Post from the coronal end.
8. Build up the core pattern with a pattern resin.
9. Remove the post/core pattern from the tooth (**Fig. 24**).
10. Insert the corresponding ParaPost Temporary Post into the prepared post space. If necessary, shorten from the apical end. Fabricate a provisional crown with self-curing resin (**Fig. 25**). To preserve the integrity of the post space, apply temporary cement **ONLY** to the margins of the crown; **DO NOT PLACE** any cement into the post space or on the post.

11. Sprue, invest and cast (**Fig. 26**).
 - A. When using a ParaPost Burnout Post, cast with any dental alloy.
 - B. When casting to a ParaPost No-Ox Precious Alloy Post, use any dental alloy with a casting temperature no higher than 2500°F (1370°C). **DO NOT USE** gypsum-bonded investment materials or those containing chloride, as they will make the prefabricated post brittle. Suggested investment materials include: Ceramigold, WM80, Ceramigold 2, Cristobalite Inlay, Hi-Temp, Cristobalite Model, Hi-Temp 2, Hydrovest, Chrome Investment X20.
12. Finish and cement the cast post and core.
13. To prepare for cementation, clean and rinse the canal and dry with Hygienic Paper Points. Use the dental cement of your choice according to the manufacturer's instructions.
14. Use a Lentulo spiral to coat the walls of the post space with cement.
15. Coat the post shaft with cement and slowly insert the post to full depth, allowing excess cement to vent.

16. In order to obtain maximum retention to crown, DO NOT POLISH the core. Fabricate the final restoration. After the cement has set, remove any excess before completing the core and fabricating the final restoration (**Fig. 27**).

Indirect Only

7. Place into the prepared post space the ParaPost Impression Post corresponding to the last drill used. If necessary, adjust the length of the Impression Post so the post does not contact the impression tray.
8. Take an impression with an elastomeric impression material (i.e., President). Make sure the anti-rotational box is filled with impression material (**Fig. 28**).
NOTE: To ensure retention of the post in the impression material when using Hydrocolloid, form a nail-head on the top of the Impression Post by pressing on the head with a hot instrument. DO NOT ALTER the plastic Impression Post when using silicone or rubber-based impression material.
9. Insert the corresponding ParaPost Temporary Post into the prepared post space. If necessary, shorten from the apical end. Fabricate a provisional crown

(**Fig. 29**). To preserve the integrity of the post space, apply temporary cement ONLY to the margins of the crown; DO NOT PLACE any cement into the post space or on the post.

10. Pour a model and fabricate a precision removable die using the Pindex® or Accutrac® System.
11. Place the corresponding ParaPost Burnout Post or ParaPost No-Ox Precious Alloy Post into the post space of the model (**Fig. 30**). If the No-Ox Post does not readily fit the post space in the die, carefully ream the post space with the corresponding ParaPost Drill in the Universal Hand Driver. If necessary, shorten the Burnout Post or No-Ox Post from the coronal end.
12. Complete the core pattern.
13. Carefully remove the post/core pattern (**Fig. 31**).
14. Sprue, invest and cast.
 - A. When using a plastic ParaPost Burnout Post, cast with any alloy.
 - B. When casting to a ParaPost No-Ox Precious Alloy Post, use any dental alloy with a casting temperature no higher than 2500°F (1370°C).

DO NOT USE gypsum-bonded investment materials or those containing chloride, as they will make the prefabricated post brittle. Suggested investment materials include: Ceramigold, WM80, Ceramigold 2, Cristobalite Inlay, Hi-Temp, Cristobalite Model, Hi-Temp 2, Hydrovest, Chrome Investment X-20.

15. Finish and cement the cast post core.
16. To prepare for cementation, clean and rinse the canal and dry with Hygenic Paper Points. Use the dental cement of your choice according to the manufacturer's instructions.
17. Use a Lentulo spiral to coat the walls of the post space with cement.
18. Coat the post shaft with cement and slowly insert the post to full depth, allowing excess cement to vent.
19. In order to obtain maximum retention to crown, DO NOT POLISH the casting. After the cement has set, remove any excess before completing the core and fabricate the final restoration (**Fig. 32**).

Cleaning and Sterilization

Drills: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush; thoroughly rinse and dry. Place drills, with a sterilization indicator, into an autoclave bag and seal. Sterilize in autoclave at 132°C for 12 minutes. Sterilizing bags should not touch the walls of the autoclave. Once the sterilization monitor indicates that sterilization was achieved, place autoclave bag, with the drills inside, in storage until ready for use. To decontaminate, soak in a 70% isopropyl alcohol solution for 10 minutes.

Temporary and No-Ox Posts: Prior to each use, sterilize by autoclave, following standard procedures.

Universal Hand Driver: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush. Sterilize by autoclave or chemical vapor sterilization.

Warning Notice: DO NOT BEND POSTS.

**Reference and/or literature available upon request.*

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Contents 7 ParaPost Drills (1 of each size)
 25 Posts (varied quantities of each size)
 1 Universal Hand Driver

Refills		ParaPostXP				ParaPostXH	
Size	Drills (3)	Stainless Steel (10)	Stainless Steel (25)	Titanium Alloy (10)	Titanium Alloy (25)	Titanium Alloy (10)	Titanium Alloy (25)
3	P-42-3	P-744-3	P-744-3B	P-784-3	P-784-3B	P-88-3	P-88-3B
4	P-42-4	P-744-4	P-744-4B	P-784-4	P-784-4B	P-88-4	P-88-4B
4.5	P-42-4.5	P-744-4.5	P-744-4.5B	P-784-4.5	P-784-4.5B	P-88-4.5	P-88-4.5B
5	P-42-5	P-744-5	P-744-5B	P-784-5	P-784-5B	P-88-5	P-88-5B
5.5	P-42-5.5	P-744-5.5	P-744-5.5B	P-784-5.5	P-784-5.5B	P-88-5.5	P-88-5.5B
6	P-42-6	P-744-6	P-744-6B	P-784-6	P-784-6B	P-88-6	P-88-6B
7	P-42-7	P-744-7	P-744-7B	P-784-7	P-784-7B	P-88-7	P-88-7B

ParaPostXP Casting Technique - Introductory Kit P-781

Contents

- 7 ParaPost Drills (1 of each size)
- 25 Plastic Impression Posts (varied quantities of each size)
- 25 Titanium Temporary Posts (varied quantities of each size)
- 25 Plastic Burnout Posts (varied quantities of each size)

Refills		ParaPostXP Casting						
Size	Drills (3)	Impression Post (20)	Titanium Post (20)	Burnout Post (10)	Burnout Post (25)	Burnout Post (100)	Nickel/Silv. Post (10)	No-Ox Post (1)
3	P-42-3	P-743-3	P-746-3	P-751-3	P-751-3B	P-751-3E	P-753-3	P-752-3
4	P-42-4	P-743-4	P-746-4	P-751-4	P-751-4B	P-751-4E	P-753-4	P-752-4
4.5	P-42-4.5	P-743-4.5	P-746-4.5	P-751-4.5	P-751-4.5B	P-751-4.5E	P-753-4.5	P-752-4.5
5	P-42-5	P-743-5	P-746-5	P-751-5	P-751-5B	P-751-5E	P-753-5	P-752-5
5.5	P-42-5.5	P-743-5.5	P-746-5.5B	P-751-5.5	P-751-5.5B	P-751-5.5E	P-753-5.5	P-752-5.5
6	P-42-6	P-743-6	P-746-6	P-751-6	P-751-6B	P-751-6E	P-753-6	P-752-6
7	P-42-7	P-743-7	P-746-7	P-751-7	P-751-7B	P-751-7E	P-753-7	P-752-7

ParaPost Refills

Refills		ParaPost		
Size	Drills (3)	Stainless Steel (10)	Stainless Steel (25)	Titanium Alloy (10)
3	P-42-3	P-44-3	P-44-3B	P-84-3
4	P-42-4	P-44-4	P-44-4B	P-84-4
4.5	P-42-4.5	P-44-4.5	P-44-4.5B	P-84-4.5
5	P-42-5	P-44-5	P-44-5B	P-84-5
5.5	P-42-5.5	P-44-5.5	P-44-5.5B	P-84-5.5
6	P-42-6	P-44-6	P-44-6B	P-84-6
7	P-42-7	P-44-7	P-44-7B	P-84-7

ParaPost Casting Technique

Refills		ParaPost Casting					
Size	Drills (3)	Impression Post (20)	Titanium Post (20)	Burnout Post (10)	Burnout Post (25)	Burnout Post (100)	Nickel/Silv. Post (10)
3	P-42-3	P-43-3	P-46-3	P-51-3	P-51-3B	P-51-3E	P-53-3
4	P-42-4	P-43-4	P-46-4	P-51-4	P-51-4B	P-51-4E	P-53-4
4.5	P-42-4.5	P-43-4.5	P-46-4.5	P-51-4.5	P-51-4.5B	P-51-4.5E	P-53-4.5
5	P-42-5	P-43-5	P-46-5	P-51-5	P-51-5B	P-51-5E	P-53-5
5.5	P-42-5.5	P-43-5.5	P-46-5.5	P-51-5.5	P-51-5.5B	P-51-5.5E	P-53-5.5
6	P-42-6	P-43-6	P-46-6	P-51-6	P-51-6B	P-51-6E	P-53-6
7	P-42-7	P-43-7	P-46-7	P-51-7	P-51-7B	P-51-7E	P-53-7


ParaPost PLUS








Refills		ParaPost Plus			
Size	Drills (3)	Stainless Steel (10)	Stainless Steel (25)	Titanium Alloy (10)	Titanium Alloy (25)
3	P-42-3	P-244-3	P-244-3B	P-284-3	P-284-3B
4	P-42-4	P-244-4	P-244-4B	P-284-4	P-284-4B
4.5	P-42-4.5	P-244-4.5	P-244-4.5B	P-284-4.5	P-284-4.5B
5	P-42-5	P-244-5	P-244-5B	P-284-5	P-284-5B
5.5	P-42-5.5	P-244-5.5	P-244-5.5B	P-284-5.5	P-284-5.5B
6	P-42-6	P-244-6	P-244-6B	P-284-6	P-284-6B
7	P-42-7	P-244-7	P-244-7B	P-284-7	P-284-7B






Approximate Diameter Comparisons

Fig. A

ParaPost® Drill				ISO D2 Instrument Measurement #/mm	Peeso Measurement #/mm		Gates-Glidden Measurement #/mm	
#	inches	mm	color					
					#1 0.72		#1 0.50 #2 0.70	
3	0.036	0.90	Brown	#55 0.87	#2 0.85	#3 0.90		
4	0.040	1.00	Yellow	#60 0.92	#2 0.85	#3 0.90		
4.5	0.045	1.14	Blue	#70 1.02	#3 1.02	#4 1.10		
5	0.050	1.25	Red	#90 1.22	#4 1.20	#4 1.10		
5.5	0.055	1.40	Purple	#100 1.32	#5 1.40	#5 1.30		
6	0.060	1.50	Black	#110 1.42	#6 1.50	#6 1.50		
7	0.070	1.75	Green	#140 1.72				

28.


29.

30.

31.

32.

ParaPost® X System

Coltène/Whaledent Inc.

235 Ascot Parkway
Cuyahoga Falls, OH 44223 / USA
Tel. USA & Canada 1 800 221 3046
+1 330 916 8800
Fax +1 330 916 7077
info@coltenewhaledent.com

Coltène/Whaledent AG

Feldwiesenstrasse 20
9450 Altstätten / Switzerland
Tel. +41 (0) 71 757 53 00
Fax +41 (0) 71 757 53 01
office@coltenewhaledent.ch

Authorized EU Representative:

Coltène/Whaledent GmbH + Co. KG

Raiffeisenstrasse 30
89129 Langenau / Germany
Tel. +49 (0) 7345 805 0
Fax +49 (0) 7345 805 201
productinfo@coltenewhaledent.de

Coltène/Whaledent Website

<http://www.coltenewhaledent.com>

Made in USA

©2003 P/N 86629J

CE
0344

coltène 
whaledent®

Caution: Federal law restricts this device to sale by or on the order of a dentist.