



# Polish with diamonds. Skip the paste.

## 3M™ Sof-Lex™ Diamond Polishing System

How much time and effort do you spend creating beautiful smiles? Whether you currently use a rubberised finishing and polishing system or an intraoral diamond polish, the process can be time-consuming and, even with your best effort, the gloss may not last. 3M has a simple solution for both problems, using two of our innovative technologies.

**Restore with 3M™ Filtek™ Supreme XTE Universal Restorative** Life-like aesthetics is just one reason why dentists use this nanocomposite. Thanks to 3M's true nanotechnology, it is easy to polish and offers excellent polish retention.<sup>1</sup>

### Polish with the 3M™ Sof-Lex™ Diamond Polishing System

Forget the messy paste. Our pre-polishing spiral prepares the restoration for final gloss, while our diamond-impregnated polishing spiral gives your restorations that glorious paste-like gloss.<sup>1</sup> The system offers the convenience of a rubberised system while also adapting to all tooth surfaces.<sup>2</sup>

You'll be happy to know that while the spirals are effective, they're also kinder to gingival tissues compared to other finishing and polishing tools<sup>3</sup> and maintain the integrity and anatomy of your restorations!

When patients leave your surgery smiling, you'll marvel at how simple it's become to create beautiful, natural-looking aesthetics.

### A difference that you can see!



**Filtek Supreme XTE universal restorative polished with Sof-Lex diamond polishing system (left)** vs. TPH Spectra® Universal Composite polished with Enhance® Finishing System and PoGo® Polishing System (right). Prepared using the manufacturer's guidelines<sup>7</sup>. Notice a clearer reflection with the Sof-Lex diamond polishing system.

You can create a diamond paste-like gloss with just two steps.



**1**  
Sof-Lex pre-polishing spiral (beige)



**2**  
Sof-Lex diamond polishing spiral (pink)

## Summary of advantages

- Imparts paste-like gloss in the convenience of a rubberised system<sup>1</sup>
- Unique, flexible shape adapts to all tooth surfaces<sup>2</sup>
- Fast and easy to use<sup>4</sup>
- Multi-use, can be sterilized and reused<sup>8</sup>
- High, long-lasting gloss when used with Filtek Supreme XTE universal restorative<sup>5</sup>

# Don't you love it when you find a product that works great, but is also cost effective?

Then you'll love this! The multipurpose 3M™ Sof-Lex™ spirals can be sterilised and reused, making them an economical choice for your practice.<sup>6</sup>

Abrasive particles are embedded throughout the rubberised spirals, therefore all surfaces of the spiral design can be used: side, top, bottom and edge.

Forget the need to stop and switch between traditional points, cups, discs and brushes, too. The flexible spirals adapt to all tooth surfaces—anterior and posterior, convex and concave.<sup>2</sup>



### 3M™ Sof-Lex™ Pre-Polishing Spiral

The pre-polishing spiral (beige) smooths and removes scratches in restorations that develop during contouring and prepares the surface for high-gloss polishing.



### 3M™ Sof-Lex™ Diamond Polishing Spiral

The diamond polishing spiral (pink) contains diamond particles to create a smooth, high gloss polish. **You're going to love how easy it is to create a beautiful smile!**

## Ordering Information

Item #	Product Information
5092-I	Sof-Lex™ Diamond Polishing System—Introductory Kit Kit includes: 5 Sof-Lex™ Pre-Polishing Spirals (Beige); 5 Sof-Lex™ Diamond Polishing Spirals (Pink); Instructions for Use; Technique Card
5090	Sof-Lex™ Pre-Polishing Spirals—Refill Includes: 15 Sof-Lex Pre-Polishing Spirals (Beige)
5091	Sof-Lex™ Diamond Polishing Spirals—Refill Includes: 15 Sof-Lex Diamond Polishing Spirals (Pink)

### References

- 1 3M Oral Care internal data. Imparts a paste-like gloss. Claim number 06549, 2015
- 2 3M Oral Care internal data. Adapts to all surfaces. Claim number 06524, 2015
- 3 3M Oral Care internal data. Kinder to Gingival Tissues than other finishing and polishing tools. Claim number 06521, 2015
- 4 3M Oral Care internal data. Efficiently achieves highly aesthetic results. Claim number 06521, 2015
- 5 3M Oral Care internal data. Gloss on composites. Claim number 06577, 2015
- 6 3M Oral Care internal data. All surfaces of the spiral can be used. Claim number 06528, 2015
- 7 3M Oral Care internal data. Marketing photo gloss comparison. Claim number 06548, 2015
- 8 3M Oral Care internal data. Reusable. Claim number 06485, 2015

**3M.com/SofLex**



3M Oral Care  
3M United Kingdom PLC  
Charnwood Campus  
10 Bakewell Road  
Loughborough  
Leicestershire  
LE11 5RB  
0845 602 5094

3M Ireland Ltd  
The Iveagh Building  
The Park  
Carrickmines  
Dublin 18  
00 353 (01) 280 3555

J322331  
GD501

TPH Spectra, Enhance and PoGo are trademarks of Dentsply International.

3M, Filtek and Sof-Lex are trademarks of the 3M Company. © 3M 2018. All rights reserved.